

For Research Use Only

NCOA2-Specific Polyclonal antibody

Catalog Number: 20031-1-AP



Basic Information

Catalog Number:

20031-1-AP

Size:

190 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006540

GeneID (NCBI):

10499

UNIPROT ID:

Q15596

Full Name:

nuclear receptor coactivator 2

Calculated MW:

159 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

NCOA2, also named as TIF2, belongs to the SRC/p160 nuclear receptor coactivator family. NCOA2 is a transcriptional coactivator for steroid receptors and nuclear receptors. NCOA2 is required with NCOA1 to control energy balance between white and brown adipose tissues. Chromosomal aberrations involving NCOA2 may be a cause of acute myeloid leukemias. This antibody is specific to NCOA2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

