

For Research Use Only

# ARMCX1 Polyclonal antibody, PBS Only

Catalog Number: 20193-1-PBS

Featured Product



## Basic Information

Catalog Number:

20193-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14120

GenBank Accession Number:

BC002691

GeneID (NCBI):

51309

UNIPROT ID:

Q9P291

Full Name:

armadillo repeat containing, X-linked

1

Calculated MW:

453 aa, 49 kDa

Observed MW:

49-55 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse ovary tissue were subjected to SDS PAGE followed by western blot with 20193-1-AP (ARMCX1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20193-1-PBS in a different storage buffer formulation.