

For Research Use Only

FAM107B Polyclonal antibody

Catalog Number:20527-1-AP



Basic Information

Catalog Number:

20527-1-AP

Size:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14389

GenBank Accession Number:

BC004872

GeneID (NCBI):

83641

UNIPROT ID:

Q9H098

Full Name:

family with sequence similarity 107, member B

Calculated MW:

131 aa, 16 kDa

Observed MW:

38-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse pancreas tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20527-1-AP (FAM107B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.