

For Research Use Only

SOX8 Polyclonal antibody

Catalog Number: 20627-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

20627-1-AP

Size:

1200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14556

GenBank Accession Number:

BC031797

GeneID (NCBI):

30812

UNIPROT ID:

P57073

Full Name:

SRX (sex determining region Y)-box 8

Calculated MW:

446 aa, 47 kDa

Observed MW:

47-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

IHC 1:20-1:200

IF 1:50-1:500

Applications

Tested Applications:

WB, IF/ICC, IHC, ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : mouse brain tissue, COLO 320 cells, rat brain tissue

IHC : human brain tissue,

IF : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Wenyan Tian	33190587	Cancer Biol Ther	IHC
Daniel A Michelson	35714609	Cell	IF
Hailin Tang	30810729	Carcinogenesis	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20627-1-AP (SOX8 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human brain using 20627-1-AP (SOX8 antibody) at dilution of 1:50 (under 10x lens).



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using 20627-1-AP (SOX8 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).