

For Research Use Only

REPS2 Polyclonal antibody

Catalog Number: 20699-1-AP **1 Publications**



Basic Information

Catalog Number:

20699-1-AP

Size:

1200 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004726

GeneID (NCBI):

9185

UNIPROT ID:

Q8NFB8

Full Name:

RALBP1 associated Eps domain containing 2

Calculated MW:

72 kDa

Observed MW:

72 kDa, 55 kDa, 30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

human

Background Information

REPS2, also named as POB1, is involved in growth factor signaling through its influence on the Ral signaling pathway. The antibody is specific to REPS2.

Notable Publications

Author	Pubmed ID	Journal	Application
Chan Wen	35690292	Cell Signal	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data