

For Research Use Only

LHX3 Polyclonal antibody

Catalog Number: 20745-1-AP **1 Publications**



Basic Information

Catalog Number:

20745-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014564

GeneID (NCBI):

8022

UNIPROT ID:

Q9UBR4

Full Name:

LIM homeobox 3

Calculated MW:

44 kDa

Observed MW:

44-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, HEK-293T cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Peng Zhao	35126731	Oncol Lett	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 20745-1-AP (LHX3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.