

For Research Use Only

# MED26 Polyclonal antibody

Catalog Number: 21043-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

21043-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15381

GenBank Accession Number:

BC110430

GeneID (NCBI):

9441

UNIPROT ID:

O95402

Full Name:

mediator complex subunit 26

Calculated MW:

600 aa, 65 kDa

Observed MW:

66-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, mouse ovary tissue, HEK-293 cells, rat brain tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xiaodong Lu	27353326	Nucleic Acids Res	WB
Gregory S Lambert	38898526	Retrovirology	WB
Gregory S Lambert	38293010	bioRxiv	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21043-1-AP (MED26 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.