

For Research Use Only

EBF3 Polyclonal antibody

Catalog Number: 21148-1-AP



Basic Information

Catalog Number:

21148-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15693

GenBank Accession Number:

BC126130

GeneID (NCBI):

253738

UNIPROT ID:

Q9H4W6

Full Name:

early B-cell factor 3

Calculated MW:

596 aa, 65 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21148-1-AP (EBF3 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.