

For Research Use Only

# GIYD2 Polyclonal antibody

Catalog Number: 21158-1-AP

Featured Product

4 Publications



## Basic Information

Catalog Number:

21158-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13094

GenBank Accession Number:

BC000803

GeneID (NCBI):

79008

UNIPROT ID:

Q9BQ83

Full Name:

GIY-YIG domain containing 2

Calculated MW:

31 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse brain tissue, fetal human brain tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhe Yang	32883681	Genes Dev	WB
Stephanie Panier	31447390	Mol Cell	
Siyuan Chen	38467834	EMBO J	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 21158-1-AP (GLYD2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.