

For Research Use Only

IQCC Polyclonal antibody

Catalog Number: 21533-1-AP



Basic Information

Catalog Number:

21533-1-AP

Size:

100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16161

GenBank Accession Number:

BC098317

GeneID (NCBI):

55721

UNIPROT ID:

Q4KMZ1

Full Name:

IQ motif containing C

Calculated MW:

466 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, MDA-MB-231 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21533-1-AP (IQCC antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.