

For Research Use Only

WDR21C Polyclonal antibody

Catalog Number: 21571-1-AP **2 Publications**



Basic Information

Catalog Number: 21571-1-AP	GenBank Accession Number: BC111778	Purification Method: Antigen affinity purification
Size: 260 µg/ml	GeneID (NCBI): 138009	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8NA75	
Isotype: IgG	Full Name: WD repeat domain 21C	
Immunogen Catalog Number: AG16098	Calculated MW: 395 aa, 44 kDa	
	Observed MW: 44-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, MCF-7 cells
Cited Applications: WB, IHC	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yu-Qin Mao	32452127	Cell Prolif	IHC
Haiyu Wang	27158335	Am J Transl Res	WB, IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 21571-1-AP (WDR21C antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.