

For Research Use Only

INMT Polyclonal antibody

Catalog Number: 21578-1-AP **1 Publications**



Basic Information

Catalog Number: 21578-1-AP	GenBank Accession Number: BC106903	Purification Method: Antigen affinity purification
Size: 133 µg/ml	GeneID (NCBI): 11185	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: O95050	
Isotype: IgG	Full Name: indolethylamine N-methyltransferase	
Immunogen Catalog Number: AG16121	Calculated MW: 263 aa, 29 kDa	
	Observed MW: 29 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse pancreas tissue,
Cited Applications: WB	
Species Specificity: human, mouse	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Huishan Sun	37202783	BMC Cancer	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with 21578-1-AP (INMT antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.