

For Research Use Only

ZNF192 Polyclonal antibody

Catalog Number: 21656-1-AP



Basic Information

Catalog Number:

21656-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16381

GenBank Accession Number:

BC096283

GeneID (NCBI):

7745

UNIPROT ID:

Q15776

Full Name:

zinc finger protein 192

Calculated MW:

578 aa, 66 kDa

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HepG2 cells, HeLa cells, mouse ovary tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 21656-1-AP (ZNF192 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.