

For Research Use Only

ZNF509 Polyclonal antibody

Catalog Number: 21772-1-AP



Basic Information

Catalog Number:

21772-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16445

GenBank Accession Number:

BC109087

GeneID (NCBI):

166793

UNIPROT ID:

Q6ZSB9

Full Name:

zinc finger protein 509

Calculated MW:

765 aa, 85 kDa

Observed MW:

85 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse skeletal muscle tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 21772-1-AP (ZNF509 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.