

For Research Use Only

IL-13 Polyclonal antibody

Catalog Number: 22085-1-AP **1 Publications**



Basic Information

Catalog Number: 22085-1-AP	GenBank Accession Number: BC096139	Purification Method: Antigen affinity purification
Size: 230 µg/ml	GeneID (NCBI): 3596	
Source: Rabbit	UNIPROT ID: P35225	
Isotype: IgG	Full Name: interleukin 13	
Immunogen Catalog Number: AG17280	Calculated MW: 146 aa, 16 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IHC

Species Specificity:
human

Cited Species:
rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Chuanfeng Fang	26722390	Int J Clin Exp Pathol	WB,IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data