

For Research Use Only

HMX2 Polyclonal antibody

Catalog Number: 23160-1-AP



Basic Information

Catalog Number:

23160-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19564

GenBank Accession Number:

BC137139

GeneID (NCBI):

3167

UNIPROT ID:

A2RU54

Full Name:

H6 family homeobox 2

Calculated MW:

273 aa, 30 kDa

Observed MW:

30 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, NCCIT cells, mouse brain tissue, rat brain tissue;

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23160-1-AP (HMX2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.