

For Research Use Only

C20orf194 Polyclonal antibody

Catalog Number: 23184-1-AP



Basic Information

Catalog Number:

23184-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19714

GenBank Accession Number:

BC122529

GeneID (NCBI):

25943

UNIPROT ID:

Q5TEA3

Full Name:

chromosome 20 open reading frame

194

Calculated MW:

1177 aa, 132 kDa

Observed MW:

132 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 23184-1-AP (C20orf194 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.