

For Research Use Only

N4BP2 Polyclonal antibody, PBS Only

Catalog Number: 23892-1-PBS



Basic Information

Catalog Number: 23892-1-PBS	GenBank Accession Number: BC126466	Purification Method: Antigen Affinity purified
Size: 1 mg/ml	GeneID (NCBI): 55728	
Source: Rabbit	UNIPROT ID: Q86UW6	
Isotype: IgG	Full Name: NEDD4 binding protein 2	
Immunogen Catalog Number: AG20943	Calculated MW: 1770 aa, 199 kDa	
	Observed MW: 199 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

N4BP2 (Nedd4 binding protein 2) is a Bcl-3 binding protein, which contains a polynucleotide kinase domain (PNK) at the N-terminus and a Small MutS Related (Smr) domain with nicking endonuclease activity near C-terminus. N4BP2 was at higher levels in a majority of tumor tissues examined, relative to paired normal tissues (PMID: 17626640).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23892-1-AP (N4BP2 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23892-1-PBS in a different storage buffer formulation.