

For Research Use Only

PTPN6 Polyclonal antibody, PBS Only

Catalog Number: 24143-1-PBS



Basic Information

Catalog Number:

24143-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21415

GenBank Accession Number:

BC002523

GeneID (NCBI):

5777

UNIPROT ID:

P29350

Full Name:

protein tyrosine phosphatase, non-receptor type 6

Calculated MW:

597 aa, 68 kDa

Observed MW:

63 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 24143-1-AP (PTPN6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24143-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 24143-1-AP (PTPN6 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 24143-1-PBS in a different storage buffer formulation.