

For Research Use Only

DENND3 Polyclonal antibody

Catalog Number: 24414-1-AP **1 Publications**



Basic Information

Catalog Number: 24414-1-AP	GenBank Accession Number: BC133023	Purification Method: Antigen affinity purification
Size: 600 µg/ml	GeneID (NCBI): 22898	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: A2RUS2	
Isotype: IgG	Full Name: DENN/MADD domain containing 3	
Immunogen Catalog Number: AG19859	Calculated MW: 1198 aa, 135 kDa	
	Observed MW: 135-145 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse stomach tissue,
Cited Applications: WB	
Species Specificity: human, mouse	
Cited Species: human	

Background Information

DENND3 may play a role in protein transport from recycling endosomes to lysosomes regulating, for instance, the degradation of the transferrin receptor

Notable Publications

Author	Pubmed ID	Journal	Application
Wenyang Li	39104482	Heliyon	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse stomach tissue were subjected to SDS PAGE followed by western blot with 24414-1-AP (DENND3 Antibody) at dilution of 1:1000 incubated at 4 degree celsius over night.