

For Research Use Only

# FAM47B Polyclonal antibody

Catalog Number: 24467-1-AP



## Basic Information

Catalog Number:

24467-1-AP

Size:

750 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21562

GenBank Accession Number:

BC035026

GeneID (NCBI):

170062

UNIPROT ID:

Q8NA70

Full Name:

family with sequence similarity 47, member B

Calculated MW:

645 aa, 74 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

## Background Information

FAM47B(Family with sequence similarity 47 member B) belongs to the FAM47 family and is expressed in testis.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24467-1-AP (FAM47B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.