

For Research Use Only

TRIM40 Polyclonal antibody

Catalog Number: 24526-1-AP



Basic Information

Catalog Number:

24526-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21649

GenBank Accession Number:

BC060785

GeneID (NCBI):

135644

UNIPROT ID:

Q6P9F5

Full Name:

tripartite motif-containing 40

Calculated MW:

258 aa, 29 kDa

Observed MW:

29-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24526-1-AP (TRIM40 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.