

For Research Use Only

C10orf125 Polyclonal antibody

Catalog Number: 24537-1-AP



Basic Information

Catalog Number: 24537-1-AP	GenBank Accession Number: BC129819	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 282969	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: A2VDF0	
Isotype: IgG	Full Name: chromosome 10 open reading frame 125	
Immunogen Catalog Number: AG21689	Calculated MW: 154 aa, 17 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human liver cancer tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 24537-1-AP (C10orf125 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 24537-1-AP (C10orf125 Antibody) at dilution of 1:50 (under 40x lens).