

For Research Use Only

SPRR2F Polyclonal antibody

Catalog Number: 24648-1-AP **1 Publications**



Basic Information

Catalog Number: 24648-1-AP	GenBank Accession Number: BC156788	Purification Method: Antigen affinity purification
Size: 1200 µg/ml	GeneID (NCBI): 6705	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q96RM1	
Isotype: IgG	Full Name: small proline-rich protein 2F	
Immunogen Catalog Number: AG20323	Calculated MW: 72 aa, 8 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : mouse ovary tissue,
Cited Applications: IF	
Species Specificity: human, mouse	
Cited Species: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xiaomu Cheng	33953164	Cell Death Dis	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse ovary tissue slide using 24648-1-AP (SPRR2F antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).