

For Research Use Only

ZNF674 Polyclonal antibody

Catalog Number: 24747-1-AP



Basic Information

Catalog Number:

24747-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21576

GenBank Accession Number:

BC157038

GeneID (NCBI):

641339

UNIPROT ID:

Q2M3X9

Full Name:

zinc finger family member 674

Calculated MW:

581 aa, 67 kDa

Observed MW:

63-67 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24747-1-AP (ZNF674 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.