

For Research Use Only

IDNK Polyclonal antibody

Catalog Number: 24831-1-AP



Basic Information

Catalog Number:

24831-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20088

GenBank Accession Number:

BC142991

GeneID (NCBI):

414328

UNIPROT ID:

Q5T6J7

Full Name:

chromosome 9 open reading frame

103

Calculated MW:

230 aa, 25 kDa

Observed MW:

26-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, HEK-293 cells, mouse kidney tissue, human testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 24831-1-AP (IDNK Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.