

For Research Use Only

ZNF195 Polyclonal antibody

Catalog Number: 24861-1-AP



Basic Information

Catalog Number:

24861-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21735

GenBank Accession Number:

BC121029

GeneID (NCBI):

7748

UNIPROT ID:

O14628

Full Name:

zinc finger protein 195

Calculated MW:

629 aa, 72 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : fetal human brain tissue, human heart tissue,
mouse skeletal muscle tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 24861-1-AP (ZNF195 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.