

For Research Use Only

FAM46C Polyclonal antibody

Catalog Number: 25038-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

25038-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19147

GenBank Accession Number:

BC131726

GeneID (NCBI):

54855

UNIPROT ID:

Q5VWP2

Full Name:

family with sequence similarity 46, member C

Calculated MW:

391 aa, 45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuan Xiao Zhu	28619709	Cancer Res	WB
Jo Kanasugi	32176823	Cancer Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data