

For Research Use Only

UHRF1BP1L Polyclonal antibody

Catalog Number: 25121-1-AP



Basic Information

Catalog Number:

25121-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19141

GenBank Accession Number:

BC127017

GeneID (NCBI):

23074

UNIPROT ID:

A0JNW5

Full Name:

UHRF1 binding protein 1-like

Calculated MW:

1464 aa, 164 kDa

Observed MW:

160-165 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue,

IP : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25121-1-AP (UHRF1BP1L Antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



IP result of anti-UHRF1BP1L (IP:25121-1-AP, 5ug; Detection:25121-1-AP 1:1000) with mouse brain tissue lysate 5000ug.