

For Research Use Only

DACH2 Polyclonal antibody

Catalog Number: 25219-1-AP



Basic Information

Catalog Number: 25219-1-AP	GenBank Accession Number: BC114950	Purification Method: Antigen affinity purification
Size: 600 µg/ml	GeneID (NCBI): 117154	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q96NX9	
Isotype: IgG	Full Name: dachshund homolog 2 (Drosophila)	
Immunogen Catalog Number: AG18403	Calculated MW: 599 aa, 65 kDa	
	Observed MW: 65 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue, HeLa cells, mouse eye tissue, mouse kidney tissue, NIH/3T3 cells
Species Specificity: human, mouse	

Background Information

Dachshund homolog 2 (DACH2), is a 599 amino acid protein, which belongs to the DACH family. DACH2 localizes in the nucleus and as a transcription factor is involved in regulation of organogenesis.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 25219-1-AP (DACH2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.