

For Research Use Only

TRIM56 Polyclonal antibody

Catalog Number: 25509-1-AP **1 Publications**



Basic Information

Catalog Number: 25509-1-AP	GenBank Accession Number: BC048194	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 81844	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q9BRZ2	
Isotype: IgG	Full Name: tripartite motif-containing 56	
Immunogen Catalog Number: AG21992	Calculated MW: 755 aa, 81 kDa	
	Observed MW: 81 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HeLa cells, Jurkat cells
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Benjamin M Liu	38556084	J Biol Chem	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25509-1-AP (TRIM56 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.