

For Research Use Only

TTC17 Polyclonal antibody

Catalog Number: 25660-1-AP



Basic Information

| | | |
|--------------------------------------|--|---|
| Catalog Number: 25660-1-AP | GenBank Accession Number: BC017238 | Purification Method: Antigen affinity purification |
| Size: 500 µg/ml | GeneID (NCBI): 55761 | Recommended Dilutions: IHC 1:50-1:500 |
| Source: Rabbit | UNIPROT ID: Q96AE7 | |
| Isotype: IgG | Full Name: tetra(ricopeptide repeat domain 17 | |
| Immunogen Catalog Number: AG22572 | Calculated MW: 1141 aa, 130 kDa | |

Applications

| | |
|------------------------------------|--|
| Tested Applications: IHC, ELISA | Positive Controls: IHC : human kidney tissue, |
| Species Specificity: human | |

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 25660-1-AP (TTC17 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 25660-1-AP (TTC17 Antibody) at dilution of 1:200 (under 40x lens).