

For Research Use Only

# SPATC1 Polyclonal antibody

Catalog Number: 25861-1-AP **2 Publications**



## Basic Information

Catalog Number:

25861-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12157

GenBank Accession Number:

BC053547

GeneID (NCBI):

375686

UNIPROT ID:

Q76KD6

Full Name:

spermatogenesis and centriole associated 1

Calculated MW:

500 aa, 52 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : human testis tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Brett Nixon	33451771	Mol Cell Proteomics	IF
Brett Nixon	30072580	Mol Cell Proteomics	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 25861-1-AP (SPATC1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.