

For Research Use Only

UT-B Polyclonal antibody

Catalog Number: 25962-1-AP **1 Publications**



Basic Information

Catalog Number:

25962-1-AP

Size:

3000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6769

GenBank Accession Number:

BC050539

GeneID (NCBI):

6563

UNIPROT ID:

Q13336

Full Name:

solute carrier family 14 (urea transporter), member 1 (Kidd blood group)

Calculated MW:

43 kDa

Observed MW:

48-55 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse kidney tissue, mouse skeletal muscle tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Liyi Xie	29279852	Biomed Res Int	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 25962-1-AP (UT-B antibody at dilution of 1:2000 incubated at room temperature for 1.5 hours.