

For Research Use Only

C16orf45 Polyclonal antibody

Catalog Number: 26037-1-AP



Basic Information

Catalog Number:

26037-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22584

GenBank Accession Number:

BC008967

GeneID (NCBI):

89927

UNIPROT ID:

Q96MC5

Full Name:

chromosome 16 open reading frame

45

Observed MW:

24-26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 26037-1-AP (C16orf45 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.