

For Research Use Only

ZNF470 Polyclonal antibody

Catalog Number: 26134-1-AP



Basic Information

Catalog Number:

26134-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23450

GenBank Accession Number:

BC136762

GeneID (NCBI):

388566

UNIPROT ID:

Q6ECI4

Full Name:

zinc finger protein 470

Observed MW:

73-83 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IF/ICC 1:800-1:3200

Applications

Tested Applications:

IF/ICC, WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : rat testis tissue,

IF/ICC : U2OS cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat testis tissue were subjected to SDS PAGE followed by western blot with 26134-1-AP (ZNF470 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using ZNF470 antibody (26134-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).