

For Research Use Only

OBFC1 Polyclonal antibody

Catalog Number: 26230-1-AP



Basic Information

Catalog Number:

26230-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23868

GenBank Accession Number:

BC017400

GeneID (NCBI):

79991

UNIPROT ID:

Q9H668

Full Name:

oligonucleotide/oligosaccharide-binding fold containing 1

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 26230-1-AP (OBFC1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.