

For Research Use Only

# C10orf46 Polyclonal antibody

Catalog Number: 26321-1-AP



## Basic Information

<b>Catalog Number:</b> 26321-1-AP	<b>GenBank Accession Number:</b> BC038294	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 900 µg/ml	<b>GeneID (NCBI):</b> 143384	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q86Y37	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 10 open reading frame 46	
<b>Immunogen Catalog Number:</b> AG23694	<b>Observed MW:</b> 41 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse liver tissue, LO2 cells
<b>Species Specificity:</b> human, mouse	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26321-1-AP (C10orf46 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.