

For Research Use Only

Enkurin Polyclonal antibody

Catalog Number: 26440-1-AP **1 Publications**



Basic Information

Catalog Number: 26440-1-AP	GenBank Accession Number: BC026165	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 219670	Recommended Dilutions: WB 1:1000-1:5000
Source: Rabbit	UNIPROT ID: Q8TC29	
Isotype: IgG	Full Name: chromosome 10 open reading frame 63	
Immunogen Catalog Number: AG23762	Observed MW: 25-30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis
Cited Applications: IF	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ting Song	36960713	J Cell Physiol	IF


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26440-1-AP (Enkurin antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.