

For Research Use Only

DRP2 Polyclonal antibody

Catalog Number: 26631-1-AP



Basic Information

Catalog Number:

26631-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24462

GenBank Accession Number:

BC111695

GeneID (NCBI):

1821

UNIPROT ID:

Q13474

Full Name:

dystrophin related protein 2

Calculated MW:

957 aa, 108 kDa

Observed MW:

100-110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 26631-1-AP (DRP2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.