

For Research Use Only

SH3GL1 Polyclonal antibody

Catalog Number: 27014-1-AP **2 Publications**



Basic Information

Catalog Number:

27014-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25550

GenBank Accession Number:

BC001270

GeneID (NCBI):

6455

UNIPROT ID:

Q99961

Full Name:

SH3-domain GRB2-like 1

Calculated MW:

41 kDa

Observed MW:

41-43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : rat brain tissue, mouse brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ronja Markworth	34486665	J Cell Sci	IF
Xiaoyi Mai	36259445	FASEB J	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 27014-1-AP (SH3GL1 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.