

For Research Use Only

# PAX2 Polyclonal antibody

Catalog Number: 27222-1-AP



## Basic Information

Catalog Number:	27222-1-AP	GenBank Accession Number:	BC148710	Purification Method:	Antigen affinity purification
Size:	N/A	GeneID (NCBI):	5076	Recommended Dilutions:	IHC 1:50-1:500
Source:	Rabbit	UNIPROT ID:	Q02962		
Isotype:	IgG	Full Name:	paired box 2		
Immunogen Catalog Number:	AG25929	Calculated MW:	431 aa, 46 kDa		

## Applications

Tested Applications:	IHC, ELISA	Positive Controls:	
Species Specificity:	human	IHC :	human renal cell carcinoma tissue,

**Note-IHC:** suggested antigen retrieval with **TE buffer pH 9.0; (\*)** Alternatively, antigen retrieval may be performed with **citrate buffer pH 6.0**

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

