

For Research Use Only

MFSD6L Polyclonal antibody

Catalog Number: 27730-1-AP



Basic Information

Catalog Number:

27730-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26876

GenBank Accession Number:

BC040487

GeneID (NCBI):

162387

UNIPROT ID:

Q8IWD5

Full Name:

major facilitator superfamily domain containing 6-like

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB: A549 cells, NCI-H1299 cells, TT cells

Background Information

MFSD6L (Major facilitator superfamily domain-containing protein 6-like) may be an integral membrane component.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27730-1-AP (MFSD6L antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.