

For Research Use Only

ASAH2 Polyclonal antibody

Catalog Number: 27742-1-AP

1 Publications



Basic Information

Catalog Number:

27742-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26908

GenBank Accession Number:

BC107105

GeneID (NCBI):

56624

UNIPROT ID:

Q9NR71

Full Name:

N-acylsphingosine amidohydrolase
(non-lysosomal ceramidase) 2

Calculated MW:

780 aa, 86 kDa

Observed MW:

100-115 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse small intestine tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Chi-Lee C Ho	38669694	Cancer Prev Res (Phila)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 27742-1-AP (ASA2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.