

For Research Use Only

FAM71F1 Polyclonal antibody

Catalog Number: 27984-1-AP



Basic Information

Catalog Number:

27984-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16525

GenBank Accession Number:

BC037248

GeneID (NCBI):

84691

UNIPROT ID:

Q96KD3

Full Name:

family with sequence similarity 71, member F1

Calculated MW:

342 aa, 39 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : A549 cells,

Background Information

FAM71F1, also named GARIN1B and FAM137A, belongs to the GARIN family. FAM71F1 is a testis-enriched protein containing a RAB2B-binding domain, a small GTPase involved in vesicle transport and membrane trafficking (PMID: 34714330).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 27984-1-AP (FAM71F1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.