

For Research Use Only

PABPN1 Polyclonal antibody

Catalog Number:28002-1-AP



Basic Information

Catalog Number:

28002-1-AP

Size:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27814

GenBank Accession Number:

BC010939

GeneID (NCBI):

8106

UNIPROT ID:

Q86U42

Full Name:

poly(A) binding protein, nuclear 1

Calculated MW:

33 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:20000-1:100000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, SH-SY5Y cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28002-1-AP (PABPN1 antibody) at dilution of 1:80000 incubated at room temperature for 1.5 hours.