

For Research Use Only

FCHSD2 Polyclonal antibody

Catalog Number: 28036-1-AP



Basic Information

Catalog Number:

28036-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27818

GenBank Accession Number:

BC010394

GeneID (NCBI):

9873

UNIPROT ID:

O94868

Full Name:

FCH and double SH3 domains 2

Observed MW:

85-95 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:12000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28036-1-AP (FCHSD2 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.