

For Research Use Only

# TMEM40 Polyclonal antibody

Catalog Number: 28066-1-AP



## Basic Information

|   |   |  |
|---|---|--|
| <b>Catalog Number:</b><br>28066-1-AP        | <b>GenBank Accession Number:</b><br>BC020658  | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>450 ug/ml                   | <b>GeneID (NCBI):</b><br>55287                | <b>Recommended Dilutions:</b><br>WB 1:1000-1:8000            |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q8WWA1                  |  |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>transmembrane protein 40 |  |
| <b>Immunogen Catalog Number:</b><br>AG27882 | <b>Observed MW:</b><br>25-30 kDa              |  |

## Applications

|  |  |
|--|--|
| <b>Tested Applications:</b><br>WB, ELISA | <b>Positive Controls:</b><br>WB : human placenta tissue, |
| <b>Species Specificity:</b><br>human     |  |

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Human placenta were subjected to SDS PAGE followed by western blot with 28066-1-AP (TMEM40 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.