

For Research Use Only

# RUNX2 Polyclonal antibody

Catalog Number: 28107-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 28107-1-AP	<b>GenBank Accession Number:</b> NM_001015051	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 270 µg/ml	<b>GeneID (NCBI):</b> 860	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q13950	
<b>Isotype:</b> IgG	<b>Full Name:</b> runt-related transcription factor 2	
<b>Immunogen Catalog Number:</b> AG27975	<b>Calculated MW:</b> 57 aa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
IF

**Species Specificity:**  
Human

**Cited Species:**  
mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Beibei Liu	37614629	Front Bioeng Biotechnol	IF

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data