

For Research Use Only

FC tag Polyclonal antibody

Catalog Number: 28108-1-AP



Basic Information

Catalog Number:

28108-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27844

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

26 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : Recombinant protein,

Background Information

This antibody detects the Fc fragment of mouse IgG2a.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28108-1-AP (FC tag antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.